

Transgelin-2

Catalog No: #AP88427



Package Size: #AP88427-1 50ug #AP88427-2 100ug #AP88427-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transgelin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Epididymis tissue protein Li 7e, SM22-alpha homolog, KIAA0120, CDABP0035
Accession No.	Uniprot ID:P37802
Uniprot	P37802
GeneID	8407;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TAGLN2

Note: This product is for in vitro research use only