

Transgelin-3

Catalog No: #AP88428



Package Size: #AP88428-1 50ug #AP88428-2 100ug #AP88428-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transgelin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Neuronal protein 22,NP22,Neuronal protein NP25,NP25
Accession No.	Uniprot ID:Q9UI15
Uniprot	Q9UI15
GeneID	29114;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TAGLN3

Note: This product is for in vitro research use only