

# Troponin T, fast skeletal muscle

Catalog No: #AP88432



Package Size: #AP88432-1 50ug #AP88432-2 100ug #AP88432-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Troponin T, fast skeletal muscle
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TnTf,Beta-TnTF,Fast skeletal muscle troponin T,fTnT
Accession No.	Uniprot ID:P45378
Uniprot	P45378
GeneID	7140;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TNNT3

Note: This product is for in vitro research use only