

## Ubiquitin-conjugating enzyme E2 A

Catalog No: #AP88443



Package Size: #AP88443-1 50ug #AP88443-2 100ug #AP88443-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ubiquitin-conjugating enzyme E2 A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	E2 ubiquitin-conjugating enzyme A, RAD6 homolog A, HR6A, Ubiquitin carrier protein A, Ubiquitin-protein ligase A, RAD6A
Accession No.	Uniprot ID:P49459
Uniprot	P49459
GeneID	7319;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: UBE2A

Note: This product is for in vitro research use only