

## Tumor necrosis factor receptor superfamily member 9

Catalog No: #AP88453



Package Size: #AP88453-1 50ug #AP88453-2 100ug #AP88453-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor receptor superfamily member 9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	4-1BB ligand receptor,CDw137,T-cell antigen 4-1BB homolog,T-cell antigen ILA,CD137,CD137,ILA
Accession No.	Uniprot ID:Q07011
Uniprot	Q07011
GeneID	3604;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TNFRSF9

Note: This product is for in vitro research use only