

Metalloproteinase inhibitor 4

Catalog No: #AP88454



Package Size: #AP88454-1 50ug #AP88454-2 100ug #AP88454-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metalloproteinase inhibitor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Tissue inhibitor of metalloproteinases 4,TIMP-4
Accession No.	Uniprot ID:O97563
Uniprot	O97563
GeneID	317694;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TIMP4

Note: This product is for in vitro research use only