

T-box transcription factor TBX21

Catalog No: #AP88455



Package Size: #AP88455-1 50ug #AP88455-2 100ug #AP88455-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T-box transcription factor TBX21
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	T-box protein 21, T-cell-specific T-box transcription factor T-bet, Transcription factor TBLYM, TBET, TBLYM
Accession No.	Uniprot ID:Q9UL17
Uniprot	Q9UL17
GeneID	30009;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TBX21

Note: This product is for in vitro research use only