

## Metalloproteinase inhibitor 3

Catalog No: #AP88457



Package Size: #AP88457-1 50ug #AP88457-2 100ug #AP88457-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Metalloproteinase inhibitor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein MIG-5, Tissue inhibitor of metalloproteinases 3, TIMP-3
Accession No.	Uniprot ID:P35625
Uniprot	P35625
GeneID	7078;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: TIMP3

Note: This product is for in vitro research use only