

Toll-like receptor 1

Catalog No: #AP88458



Package Size: #AP88458-1 50ug #AP88458-2 100ug #AP88458-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Toll-like receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Toll/interleukin-1 receptor-like protein,TIL,CD281
Accession No.	Uniprot ID:Q9EPQ1
Uniprot	Q9EPQ1
GeneID	21897;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tlr1

Note: This product is for in vitro research use only