

E3 ubiquitin-protein ligase SMURF2

Catalog No: #AP88464



Package Size: #AP88464-1 50ug #AP88464-2 100ug #AP88464-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase SMURF2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HECT-type E3 ubiquitin transferase SMURF2, SMAD ubiquitination regulatory factor 2, SMAD-specific E3 ubiquitin-protein ligase 2
Accession No.	Uniprot ID:A2A5Z6
Uniprot	A2A5Z6
GeneID	66313;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Smurf2

Note: This product is for in vitro research use only