

Tetraspanin-14

Catalog No: #AP88466



Package Size: #AP88466-1 50ug #AP88466-2 100ug #AP88466-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tetraspanin-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tspan-14, Transmembrane 4 superfamily member 14, D14Ert226e, Tm4sf14
Accession No.	Uniprot ID:Q8QZY6
Uniprot	Q8QZY6
GeneID	52588;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tspan14

Note: This product is for in vitro research use only