

IGFL2 Antibody

Catalog No: #44519

Package Size: #44519-1 50ul #44519-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IGFL2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	IGFL2 Antibody detects endogenous levels of total IGFL2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IGFL2
Target Name	IGFL2
Other Names	IGFL2, VPRI645, IGF-like family member 2, UNQ645
Accession No.	Swiss-Prot#: Q6UWQ7NCBI Gene ID:147920
Uniprot	Q6UWQ7
GeneID	147920;
Calculated MW	13 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only