

Transcription factor E2-alpha

Catalog No: #AP88470

Package Size: #AP88470-1 50ug #AP88470-2 100ug #AP88470-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor E2-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Immunoglobulin enhancer-binding factor E12/E47, Pancreas specific transcription factor 1c, Transcription factor 3, TCF-3, Transcription regulator Pan, Pan, Ptf1c, Tcfe2a
Accession No.	Uniprot ID:P21677
Uniprot	P21677
GeneID	171046;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tcf3

Note: This product is for in vitro research use only