

Sushi repeat-containing protein SRPX

Catalog No: #AP88472



Package Size: #AP88472-1 50ug #AP88472-2 100ug #AP88472-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sushi repeat-containing protein SRPX
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Down-regulated by v-SRC,DRS,Drs
Accession No.	Uniprot ID:Q63769
Uniprot	Q63769
GeneID	64316;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SrpX

Note: This product is for in vitro research use only