

Serine/arginine-rich splicing factor 1

Catalog No: #AP88473



Package Size: #AP88473-1 50ug #AP88473-2 100ug #AP88473-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/arginine-rich splicing factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Splicing factor, arginine/serine-rich 1,SFRS1
Accession No.	Uniprot ID:Q0VCY7
Uniprot	Q0VCY7
GeneID	615796;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SRSF1

Note: This product is for in vitro research use only