

# Tumor necrosis factor receptor superfamily member 1A

Catalog No: #AP88475

Package Size: #AP88475-1 50ug #AP88475-2 100ug #AP88475-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tumor necrosis factor receptor superfamily member 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tumor necrosis factor receptor 1,TNF-R1,Tumor necrosis factor receptor type I,TNF-RI,p55,p60,CD120a,Tnfr-1,Tnfr1
Accession No.	Uniprot ID:P25118
Uniprot	P25118
GeneID	21937;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tnfrsf1a

Note: This product is for in vitro research use only