

Tumor necrosis factor receptor superfamily member 6

Catalog No: #AP88476



Package Size: #AP88476-1 50ug #AP88476-2 100ug #AP88476-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Apo-1 antigen, Apoptosis-mediating surface antigen FAS, FASLG receptor, CD95, APT1, TNFRSF6
Accession No.	Uniprot ID:P51867
Uniprot	P51867
GeneID	282488;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FAS

Note: This product is for in vitro research use only