

Thrombopoietin

Catalog No: #AP88477



Package Size: #AP88477-1 50ug #AP88477-2 100ug #AP88477-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thrombopoietin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C-mpl ligand,ML,Megakaryocyte colony-stimulating factor,Megakaryocyte growth and development factor,MGDF,Myeloproliferative leukemia virus oncogene ligand,MGDF
Accession No.	Uniprot ID:P40225
Uniprot	P40225
GeneID	7066;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:THPO

Note: This product is for in vitro research use only