

Receptor-interacting serine/threonine-protein kinase 1

Catalog No: #AP88478



Package Size: #AP88478-1 50ug #AP88478-2 100ug #AP88478-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Receptor-interacting serine/threonine-protein kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cell death protein RIP, Receptor-interacting protein 1, RIP-1, Serine/threonine-protein kinase RIP, RIP, RIP1
Accession No.	Uniprot ID:Q13546
Uniprot	Q13546
GeneID	8737;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RIPK1

Note: This product is for in vitro research use only