

Apoptosis regulatory protein Siva

Catalog No: #AP88479



Package Size: #AP88479-1 50ug #AP88479-2 100ug #AP88479-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apoptosis regulatory protein Siva
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CD27-binding protein,CD27BP,Siva
Accession No.	Uniprot ID:P59692
Uniprot	P59692
GeneID	362791;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Siva1

Note: This product is for in vitro research use only