

IGFL4 Antibody

Catalog No: #44520

Package Size: #44520-1 50ul #44520-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IGFL4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	IGFL4 Antibody detects endogenous levels of total IGFL4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IGFL4
Target Name	IGFL4
Other Names	Insulin growth factor-like family member 4
Accession No.	Swiss-Prot#: Q6B9Z1NCBI Gene ID:444882
Uniprot	Q6B9Z1
GeneID	444882;
Calculated MW	14 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only