

## Structural maintenance of chromosomes protein 4

Catalog No: #AP88480



Package Size: #AP88480-1 50ug #AP88480-2 100ug #AP88480-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Structural maintenance of chromosomes protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SMC protein 4,Chromosome-associated polypeptide C,hCAP-C,XCAP-C homolog,CAPC,SMC4L1
Accession No.	Uniprot ID:Q9NTJ3
Uniprot	Q9NTJ3
GeneID	10051;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SMC4

Note: This product is for in vitro research use only