

TYRO protein tyrosine kinase-binding protein

Catalog No: #AP88482



Package Size: #AP88482-1 50ug #AP88482-2 100ug #AP88482-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TYRO protein tyrosine kinase-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DNAX-activation protein 12,Killer-activating receptor-associated protein,KAR-associated protein,Dap12,Karap
Accession No.	Uniprot ID:O54885
Uniprot	O54885
GeneID	22177;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tyrobp

Note: This product is for in vitro research use only