

Thymidylate synthase

Catalog No: #AP88483



Package Size: #AP88483-1 50ug #AP88483-2 100ug #AP88483-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thymidylate synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TS,TS,OK/SW-cl.29
Accession No.	Uniprot ID:P04818
Uniprot	P04818
GeneID	7298;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TYMS

Note: This product is for in vitro research use only