

Hepatitis A virus cellular receptor 2 homolog

Catalog No: #AP88486



Package Size: #AP88486-1 50ug #AP88486-2 100ug #AP88486-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepatitis A virus cellular receptor 2 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HAVcr-2,T-cell immunoglobulin and mucin domain-containing protein 3,TIMD-3,T-cell membrane protein 3,TIM-3,Tim3,Timd3,TIM-3
Accession No.	Uniprot ID:P0C0K5
Uniprot	P0C0K5
GeneID	363578;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Havcr2

Note: This product is for in vitro research use only