

Tumor necrosis factor

Catalog No: #AP88489



Package Size: #AP88489-1 50ug #AP88489-2 100ug #AP88489-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Cachectin, TNF-alpha, Tumor necrosis factor ligand superfamily member 2, TNF-a, TNFA, TNFSF2
Accession No.	Uniprot ID:P51742
Uniprot	P51742
GeneID	403922;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TNF

Note: This product is for in vitro research use only