

SUMO-conjugating enzyme UBC9

Catalog No: #AP88498



Package Size: #AP88498-1 50ug #AP88498-2 100ug #AP88498-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SUMO-conjugating enzyme UBC9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	RING-type E3 SUMO transferase UBC9,SUMO-protein ligase,Ubiquitin carrier protein 9,mUBC9,Ubiquitin carrier protein 1,Ubiquitin-conjugating enzyme E2 I,Ubiquitin-protein ligase I,Ubc9,Ubce2i,Ubce9
Accession No.	Uniprot ID:P63280
Uniprot	P63280
GeneID	22196;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ube2i

Note: This product is for in vitro research use only