

Protransforming growth factor alpha

Catalog No: #AP88501



Package Size: #AP88501-1 50ug #AP88501-2 100ug #AP88501-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protransforming growth factor alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Accession No.	Uniprot ID:Q06922
Uniprot	Q06922
GeneID	397484;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TGFA

Note: This product is for in vitro research use only