

Transforming growth factor beta-2

Catalog No: #AP88502



Package Size: #AP88502-1 50ug #AP88502-2 100ug #AP88502-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transforming growth factor beta-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TGF-beta-2,Milk growth factor,MGF
Accession No.	Uniprot ID:P21214
Uniprot	P21214
GeneID	534069;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TGFB2

Note: This product is for in vitro research use only