

# Protein SSX2

Catalog No: #AP88506



Package Size: #AP88506-1 50ug #AP88506-2 100ug #AP88506-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein SSX2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cancer/testis antigen 5.2,CT5.2,Synovial sarcoma, X breakpoint 2,Tumor antigen HOM-MEL-40,SSX2A
Accession No.	Uniprot ID:Q16385
Uniprot	Q16385
GeneID	6757;727837;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SSX2

Note: This product is for in vitro research use only