

Serglycin

Catalog No: #AP88509



Package Size: #AP88509-1 50ug #AP88509-2 100ug #AP88509-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serglycin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mastocytoma proteoglycan core protein,Secretory granule proteoglycan core protein,gp600,Prg,Prg1
Accession No.	Uniprot ID:P13609
Uniprot	P13609
GeneID	19073;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Srgn

Note: This product is for in vitro research use only