

IL31R Antibody

Catalog No: #44523

Package Size: #44523-1 50ul #44523-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IL31R Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	IL31R Antibody detects endogenous levels of total IL31R
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IL31R
Target Name	IL31R
Other Names	IL31RA, CRL, CRL3, Cytokine receptor-like 3, GLM-R, GPL, IL-31R subunit alpha, IL-31R-alpha, IL-31RA, Interleukin 31 receptor A, GLMR, HGLM-R, IL-31 receptor subunit alpha, PLCA2, ZcytoR17, Class I cytokine receptor, gp130-like monocyte receptor, gp1 ...
Accession No.	Swiss-Prot#: Q8NI17NCBI Gene ID:133396
Uniprot	Q8NI17
GeneID	133396;
Calculated MW	83 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of extracts from Jurkat cells, using IL-31R Polyclonal Antibody. Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only