

Serine/arginine-rich splicing factor 2

Catalog No: #AP88511



Package Size: #AP88511-1 50ug #AP88511-2 100ug #AP88511-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/arginine-rich splicing factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein PR264, Splicing component, 35 kDa, Splicing factor SC35, SC-35, Splicing factor, arginine/serine-rich 2, SFRS2
Accession No.	Uniprot ID:Q01130
Uniprot	Q01130
GeneID	6427;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SRSF2

Note: This product is for in vitro research use only