

Metalloproteinase inhibitor 2

Catalog No: #AP88512



Package Size: #AP88512-1 50ug #AP88512-2 100ug #AP88512-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metalloproteinase inhibitor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Collagenase inhibitor, Tissue inhibitor of metalloproteinases 2, TIMP-2
Accession No.	Uniprot ID:P16368
Uniprot	P16368
GeneID	282093;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TIMP2

Note: This product is for in vitro research use only