

Toll/interleukin-1 receptor domain-containing adapter protein

Catalog No: #AP88513

Package Size: #AP88513-1 50ug #AP88513-2 100ug #AP88513-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Toll/interleukin-1 receptor domain-containing adapter protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TIR domain-containing adapter protein,Adaptor protein Wyatt,MyD88 adapter-like protein,MyD88-2,MAL
Accession No.	Uniprot ID:P58753
Uniprot	P58753
GeneID	114609;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TIRAP

Note: This product is for in vitro research use only