

Microtubule-associated protein tau

Catalog No: #AP88518



Package Size: #AP88518-1 50ug #AP88518-2 100ug #AP88518-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Microtubule-associated protein tau
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Neurofibrillary tangle protein, Paired helical filament-tau, PHF-tau, Mtap, Tau
Accession No.	Uniprot ID:P10637
Uniprot	P10637
GeneID	17762;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mapt

Note: This product is for in vitro research use only