

Vascular endothelial growth factor D

Catalog No: #AP88524



Package Size: #AP88524-1 50ug #AP88524-2 100ug #AP88524-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vascular endothelial growth factor D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	VEGF-D,c-Fos-induced growth factor,FIGF,FIGF
Accession No.	Uniprot ID:O43915
Uniprot	O43915
GeneID	2277;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:VEGFD

Note: This product is for in vitro research use only