

KRT81/86 Antibody

Catalog No: #44525



Package Size: #44525-1 50ul #44525-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KRT81/86 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	KRT81/86 Antibody detects endogenous levels of total KRT81/86
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human KRT81/86
Target Name	KRT81/86
Other Names	KRTHB6 KRTHB1,MLN137
Accession No.	Swiss-Prot#: O43790/Q14533NCBI Gene ID:3892/3887
Uniprot	O43790
GeneID	3892;
Calculated MW	54 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only