

Tandem C2 domains nuclear protein

Catalog No: #AP88532



Package Size: #AP88532-1 50ug #AP88532-2 100ug #AP88532-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tandem C2 domains nuclear protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Membrane targeting tandem C2 domain-containing protein 1, Tandem C2 protein in nucleus, Tac2-N, Mtac2d1
Accession No.	Uniprot ID:Q91XT6
Uniprot	Q91XT6
GeneID	74413;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tc2n

Note: This product is for in vitro research use only