

Tenomodulin

Catalog No: #AP88534



Package Size: #AP88534-1 50ug #AP88534-2 100ug #AP88534-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tenomodulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	TeM, Chondromodulin-1-like protein, ChM1L, Myodulin, Chm1l
Accession No.	Uniprot ID:Q9ESC2
Uniprot	Q9ESC2
GeneID	64104;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tnmd

Note: This product is for in vitro research use only