

Tumor necrosis factor receptor superfamily member 6

Catalog No: #AP88535



Package Size: #AP88535-1 50ug #AP88535-2 100ug #AP88535-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Apo-1 antigen,Apoptosis-mediating surface antigen FAS,FASLG receptor,CD95,APT1,FAS1,TNFRSF6
Accession No.	Uniprot ID:P25445
Uniprot	P25445
GeneID	355;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FAS

Note: This product is for in vitro research use only