

Src-like-adapter 2

Catalog No: #AP88536



Package Size: #AP88536-1 50ug #AP88536-2 100ug #AP88536-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Src-like-adapter 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Src-like adapter protein 2,SLAP-2
Accession No.	Uniprot ID:Q8R4L0
Uniprot	Q8R4L0
GeneID	77799;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SlA2

Note: This product is for in vitro research use only