

Transcobalamin-1

Catalog No: #AP88542



Package Size: #AP88542-1 50ug #AP88542-2 100ug #AP88542-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcobalamin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TC-1,Haptocorrin,HC,Protein R,Transcobalamin I,TC I,TC1
Accession No.	Uniprot ID:P20061
Uniprot	P20061
GeneID	6947;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TCN1

Note: This product is for in vitro research use only