

Telomerase protein component 1

Catalog No: #AP88544



Package Size: #AP88544-1 50ug #AP88544-2 100ug #AP88544-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Telomerase protein component 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Telomerase-associated protein 1, Telomerase protein 1, p240, p80 telomerase homolog, TLP1, TP1
Accession No.	Uniprot ID:Q99973
Uniprot	Q99973
GeneID	7011;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TEP1

Note: This product is for in vitro research use only