

Lymphotoxin-alpha

Catalog No: #AP88545



Package Size: #AP88545-1 50ug #AP88545-2 100ug #AP88545-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lymphotoxin-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	LT-alpha,TNF-beta,Tumor necrosis factor ligand superfamily member 1,TNFB,TNFSF1
Accession No.	Uniprot ID:Q5WR07
Uniprot	Q5WR07
GeneID	607183;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LTA

Note: This product is for in vitro research use only