

Telomerase protein component 1

Catalog No: #AP88548



Package Size: #AP88548-1 50ug #AP88548-2 100ug #AP88548-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Telomerase protein component 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	rTLP1,Telomerase-associated protein 1,Telomerase protein 1,p230,p240,p80 telomerase homolog,Tlp1
Accession No.	Uniprot ID:O08653
Uniprot	O08653
GeneID	64523;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tep1

Note: This product is for in vitro research use only