

Serine/arginine-rich splicing factor 2

Catalog No: #AP88550



Package Size: #AP88550-1 50ug #AP88550-2 100ug #AP88550-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/arginine-rich splicing factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Splicing component, 35 kDa, Splicing factor SC35, SC-35, Splicing factor, arginine/serine-rich 2, SFRS2
Accession No.	Uniprot ID:Q06A98
Uniprot	Q06A98
GeneID	768117;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SRSF2

Note: This product is for in vitro research use only