

Serine/threonine-protein kinase STK11

Catalog No: #AP88551



Package Size: #AP88551-1 50ug #AP88551-2 100ug #AP88551-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase STK11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Liver kinase B1 homolog,LKB1,Lkb1
Accession No.	Uniprot ID:Q9WTK7
Uniprot	Q9WTK7
GeneID	20869;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Stk11

Note: This product is for in vitro research use only