

Tachykinin-3

Catalog No: #AP88554



Package Size: #AP88554-1 50ug #AP88554-2 100ug #AP88554-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tachykinin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Neurokinin B-like protein Zneurok1,Preprotachykinin-B,PPT-B,Nknb,Tac2
Accession No.	Uniprot ID:P55099
Uniprot	P55099
GeneID	21334;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tac3

Note: This product is for in vitro research use only