

Thrombomodulin

Catalog No: #AP88555



Package Size: #AP88555-1 50ug #AP88555-2 100ug #AP88555-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thrombomodulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	TM,CD141
Accession No.	Uniprot ID:Q5W7P8
Uniprot	Q5W7P8
GeneID	474355;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:THBD

Note: This product is for in vitro research use only